CLAIMS

ASSAY FOR ANTI TRANSGLUTAMINASE ANTIBODIES DETECTION USEFUL IN CELICAC DISEASE DIAGNOSIS

- 1. An assay to detect IgA or IgG anti-transglutaminase antibodies in liquid samples which comprises:
 - a) Immunocomplex formation between antibodies in the samples and the antigen tissue transglutaminase, obtained from natural sources or by recombinant DNA technology. The antigen is conjugated to a colored substance and deposited onto an inert support which allows the release of the conjugated when it comes into contact with a liquid sample.
 - b) Reaction of the immunocomplexes described in a) with the same antigen tissue transglutaminase adsorbed onto a reactive zone of a membrane promoting its deposition on this zone. We correlation step 5

The promoting its deposition on this zone. We correlation step: compore 2. An assay to detect IgA or IgG anti-transglutaminase antibodies according to claim 1, sund wherein the adsorbed antigen is fixed onto a nitrocellulose or nylon membrane with 5 to 10 µm pore size that allows a lateral flow of reactants.

3. An assay to detect IgA or IgG anti-transglutaminase antibodies according to claim 1, wherein the colored substance conjugated to the antigen is colloidal gold or colored latex particles.

- 4. In assay to detect IgA or IgG anti-transglutaminase antibodies according to slaim 1, wherein the excess of the conjugated antigen reacts with a reagent adsorbed onto another zone of the same membrane resulting in a second colored signal that can be used to check the performance of the assay.
- 7 5. Procedure according to claim 1 wherein the assay detects IgA or IgG antitransglutaminase antibodies in samples of human serum, plasma, or blood in just onestep.